Infuse™ Bone Graft

Important Medical Information

0381204E Rev. D

2019-11-04

Medtronic

ENGLISH

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician with appropriate training.

DESCRIPTION

Infuse™ Bone Graft consists of two components – a recombinant human bone morphogenetic protein solution and a carrier/scaffold for the bone morphogenetic protein solution and resulting bone. **These components must be used as a system. The bone morphogenetic protein solution component must not be used without the carrier/scaffold component or with a carrier/scaffold component different from the one described in this document.**

Infuse™ Bone Graft consists of recombinant human Bone Morphogenetic Protein-2 (rhBMP-2, known as dibotermin alfa) placed on an absorbable collagen sponge (ACS). Infuse™ Bone Graft induces new bone tissue at the site of implantation. Based on data from non-clinical studies, the bone formation process develops from the outside of the implant towards the center until the entire device is replaced by trabecular bone.

rhBMP-2 is the active agent in Infuse™ Bone Graft. rhBMP-2 is a disulfide-linked dimeric protein molecule with two major subunit species of 114 and 131 amino acids. Each subunit is glycosylated at one site with high-mannose-type glycans. rhBMP-2 is produced by a genetically engineered Chinese hamster ovary cell line.

rhBMP-2 and excipients are lyophilized. Upon reconstitution, each milliliter of rhBMP-2 solution contains: 1.5mg of rhBMP-2; 5.0 mg sucrose, NF; 25mg glycine, USP; 3.7mg L-glutamic acid, FCC; 0.1mg sodium chloride, USP; 0.1mg polysorbate 80, NF; and 1.0mL of sterile water. The reconstituted rhBMP-2 solution has a pH of 4.5, and is clear, colorless to slightly yellow and essentially free from plainly visible particulate matter.

The ACS is a soft, white, pliable, absorbent implantable matrix for rhBMP-2. ACS is made from bovine Type I collagen obtained from the deep flexor (Achilles) tendon. The ACS acts as a carrier for the rhBMP-2 and acts as a scaffold for new bone formation.

Each kit contains all the components necessary to prepare Infuse™ Bone Graft: the rhBMP-2 which must be reconstituted, sterile water, absorbable collagen sponges, syringes with needles, this package insert and instructions for preparation.

The rhBMP-2 is provided as a lyophilized powder in vials delivering 12mg of protein. After appropriate reconstitution, the concentration of rhBMP-2 is 1.5mg/mL. The solution is then applied to the provided absorbable collagen sponge. Infuse™ Bone Graft is prepared at the time of surgery and allowed a prescribed amount of time (no less than 15 minutes) before placement at the fracture site. The Instructions for Preparation contain complete details on preparation of Infuse™ Bone Graft.

INDICATIONS

Infuse™ Bone Graft is indicated for treating acute, open tibial shaft fractures that have been stabilized with IM nail fixation after appropriate wound management. Infuse™ Bone Graft must be applied within 14 days after the initial fracture. Prospective patients should be skeletally mature.

CONTRAINDICATIONS

- Infuse™ Bone Graft is contraindicated for patients with a known hypersensitivity to recombinant human Bone Morphogenetic Protein-2, bovine Type I collagen or to other components of the formulation.
- Infuse™ Bone Graft should not be used in the vicinity of a resected or extant tumor, in patients with any active malignancy or patients undergoing treatment for a malignancy.
- Infuse™ Bone Graft should not be used in patients who are skeletally immature (<18 years of age or no radiographic evidence of epiphyseal closure).
- Infuse™ Bone Graft should not be used in patients with an inadequate neurovascular status (e.g. high risk of amputation).
- Infuse™ Bone Graft should not be used in patients with compartment syndrome of the affected limb.
- Infuse™ Bone Graft should not be used in pregnant women. The potential effects of rhBMP-2 on the human fetus have not been evaluated.
- Infuse™ Bone Graft should not be implanted in patients with an active infection at the operative site.

WARNINGS

- In an experimental rabbit study, rhBMP-2 has been shown to elicit antibodies that are capable of crossing the placenta. Reduced ossification of the frontal and parietal bones of the skull was noted infrequently (<3%) in fetuses of rabbit dams immunized to rhBMP-2; however, there was no effect noted in limb bud development. There are no adequate and well-controlled studies in human pregnant women. Women of child bearing potential should be warned by their surgeon of potential risk to a fetus and informed of other possible orthopedic treatments.
- Women of childbearing potential should be advised that antibody formation to rhBMP-2 or its influence on human fetal development has not been assessed. In the clinical trial supporting the safety and effectiveness of the Infuse™ Bone Graft in tibia fractures, 9/149 (6.0%) patients treated with Infuse™ Bone Graft and 1/150 (0.7%) patients treated without exposure to rhBMP-2 developed antibodies to rhBMP-2. The effect of maternal antibodies to rhBMP-2, as might be present for several months following device implantation, on the unborn fetus is unknown. Additionally, it is unknown whether fetal expression of BMP-2 could re-expose mothers who were previously antibody positive. Theoretically, re-exposure may elicit a more powerful immune response to BMP-2 with possible adverse consequences for the fetus. However, pregnancy did not lead to an increase in antibodies in the rabbit study. Studies in genetically altered mice indicate that BMP-2 is critical to fetal development and that a lack of BMP-2 activity may cause neonatal death or birth defects. It is not known if anti-BMP-2 antibodies may affect fetal development or the extent to which these antibodies may reduce BMP-2 activity.
- Infuse™ Bone Graft should not be used immediately prior to or during pregnancy. Women of childbearing potential should be advised not to become pregnant for one year following treatment with the Infuse™ Bone Graft.
- The safety and effectiveness of the Infuse™ Bone Graft in nursing mothers has not been established. It is not known if BMP-2 is excreted in human milk.

- Infuse™ Bone Graft should not be used in patients suspected of having a malignancy at the site of application.
- The safety and effectiveness of Infuse™ Bone Graft for non-acute fractures, with forms of internal fracture fixation other than IM nails, implanted at locations other than the tibial shaft, or used in surgical techniques other than open reduction and internal fixation after appropriate wound management have not been established.
- Inappropriate use of the product, such as preparing it differently than prescribed, compressing the rhBMP-2/ACS implant more than necessary, or overfilling the volume intended for new bone formation, may change the concentration of the rhBMP-2, which may inhibit the ability of the rhBMP-2/ACS to convert to bone and/or cause complications. Such use of the rhBMP-2/ACS implant may result in radiographic evidence of resorption, fluid formation, and/or edema. These findings may be asymptomatic or symptomatic.
- The formation of fluid collections (sometimes encapsulated) in some cases resulted in pain, which may require clinical intervention (aspiration and/or surgical removal) if symptoms persist. Many of these reports have occurred when rhBMP-2/ACS was used in conjunction with unapproved approaches/devices or in a manner inconsistent with the instructions for use.
- In a clinical study in which the intramedullary canal was reamed to cortical chatter, an increased rate of infection was observed in those patients treated with rhBMP-2/ACS compared to a control group that did not receive rhBMP-2/ACS. Because of this, the use of rhBMP-2/ACS with reamed IM nails is not recommended.
- While there are currently anecdotal and literature evidence to suggest that volume overfilling and/or hyperconcentration of the rhBMP-2 solution may lead to fluid formation and/or edema, animal models for scientifically evaluating these events do not presently exist. A sheep model developed to test the hypothesis that volume overfilling and/or hyperconcentration of the rhBMP-2 solution results in radiographic evidence of bone resorption has preliminarily been evaluated and appears to be supportive of the hypothesized mechanism.
- Placement of rhBMP-2/ACS can cause initial resorption of trabecular bone that may be transient.

PRECAUTIONS

General

- The safety and effectiveness of repeat applications of Infuse™ Bone Graft has not been established.
- Long-bone fracture and soft-tissue management procedures should be based on standard practice including control of infection. Physicians should achieve mechanical stability before implanting Infuse™ Bone Graft. Infuse™ Bone Graft should not be used to fill space in the presence of compressive forces.
- Infuse™ Bone Graft should only be used by surgeons who are experienced in treating acute open tibial shaft fractures involving IM nail stabilization and have undergone adequate training with this device.
- A single package of Infuse™ Bone Graft should be used at the fracture site.
- Infuse™ Bone Graft is intended for single use only. Discard unused product and use a new device for subsequent applications. Infuse™ Bone Graft must not be sterilized by the hospital.
- Prior to use, inspect the packaging, vials and stoppers for visible damage. If damage is visible, do not use the product. Retain the packaging and vials and contact a Medtronic representative.
- Do not use after the printed expiration date on the label.

Hepatic and Renal Impairment

The safety and effectiveness of Infuse™ Bone Graft in patients with hepatic or renal impairment has not been established. Pharmacokinetic studies of rhBMP-2 indicate that the renal and hepatic systems are involved with its clearance.

Geriatrics

Clinical studies of Infuse™ Bone Graft did not include sufficient numbers of patients 65 years and older to determine whether they respond differently from younger subjects.

Bone Formation

- The safety and effectiveness of Infuse™ Bone Graft has not been demonstrated in patients with metabolic bone diseases, or in pathological fractures such as those observed in Paget's disease or in metastatic bone disease.
- The potential for ectopic, heterotopic or undesirable exuberant bone formation exists.

Antibody Formation/Allergic Reactions

- The safety and effectiveness of Infuse™ Bone Graft has not been demonstrated in patients with autoimmune disease.
- The safety and effectiveness of Infuse™ Bone Graft has not been demonstrated in patients with immunosuppressive disease or suppressed immune systems resulting from radiation therapy, chemotherapy, steroid therapy or other treatments.

Immunogenicity

- As with all therapeutic proteins, there is a potential for immune responses to be generated to Infuse™ Bone Graft. The immune response to Infuse™ Bone Graft was evaluated in 149 investigational patients and 150 control patients receiving treatment for acute open tibial shaft fractures stabilized with IM nails.
 - *Anti-rhBMP-2 antibodies:* 9/149 (6%) patients receiving the Infuse™ Bone Graft component developed antibodies vs. 1/150 (0.7%) in the control group.
 - *Anti-bovine Type I collagen antibodies:* 29/149 (20%) of patients receiving the Infuse™ Bone Graft component developed antibodies to bovine Type I collagen vs. 9/150 (6%) of control patients. No patients in either group developed anti-human Type I collagen antibodies.
 - The presence of antibodies to rhBMP-2 was not associated with immune mediated adverse events such as allergic reactions. The neutralizing capacity of antibodies to rhBMP-2 is not known.
- The incidence of antibody detection is highly dependent on the sensitivity and specificity of the assay. Additionally, the incidence of antibody detection may be influenced by several factors including sample handling, concomitant medications and underlying disease. For these reasons, comparison of the incidence of antibodies to Infuse™ Bone Graft with the incidence of antibodies to other products may be misleading.

ADVERSE EVENTS

Table 1 describes the adverse events observed in the clinical trial used to support approval of the product. Two Infuse™ Bone Graft doses, 0.75mg/ml and 1.5mg/ml, were evaluated. Infuse™ Bone Graft with IM nail stabilization was implanted in 300 investigational patients (149 in the 1.50mg/ml and 151 in the 0.75mg/ml groups) compared to IM stabilization alone in 150 control patients. Adverse event rates presented are based on the number of patients having at least one occurrence for a particular adverse event divided by the total number of patients in that treatment group.

Table 1: Adverse events reported for patients enrolled in clinical trial supporting approval

													# (%) of Patients Total Adverse Events		
	1st quarter postop			2nd quarter postop			3rd quarter postop			4th quarter postop			Ctrl (n=150) (%)	Low dose (n=151) (%)	Inv (n=149) (%)
	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv			
Number of patients	150	150	149	147	144	146	144	144	144	138	142	141	150 (100)	150 (99)	149 (100)
Abnormal healing															
Surgical Site	64	60	53	19	19	13	7	4	7	4	2	3	65 ¹ (43) ²	62 (41)	57 (38)
													94 ³	85	76
Other locations	9	13	23	0	1	0	1	0	1	0	0	2	10 (7)	10 (7)	13 (9)
													10	14	26
Abnormal lab tests															
Alkaline phosphatase increased	8	8	3	0	0	0	0	0	0	0	0	0	8 (5)	8 (5)	3 (2)
													8	8	3
Amylase increased	5	23	11	0	0	0	0	0	0	0	0	0	5 (3)	20 (13)	10 (7)
													5	23	11
Bilirubinemia	8	9	11	0	0	0	0	0	0	0	0	0	7 (5)	9 (6)	11 (7)
													8	9	11
BUN increased	0	3	1	0	0	0	0	0	0	0	0	0	0 (0)	3 (2)	1 (1)
													0	3	1
Creatinine clearance decreased	2	1	2	0	0	0	0	0	0	0	0	0	2 (1)	1 (1)	2 (1)
													2	1	2
Gamma glutamyl transpeptidase increased	0	1	2	0	0	0	0	0	0	0	0	0	0 (0)	1 (1)	2 (1)
													0	1	2
Hypercalcemia	2	1	0	0	0	0	0	0	0	0	0	0	2 (1)	1 (1)	0 (1)
													2	1	0
Hyperkalemia	2	3	3	0	0	0	0	0	0	0	0	0	2 (1)	3 (2)	3 (2)
													2	3	3

Table 1: Adverse events reported for patients enrolled in clinical trial supporting approval (continued)

													# (%) of Patients Total Adverse Events					
	1st quarter postop			2nd quarter postop			3rd quarter postop			4th quarter postop								
	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl (n=150) (%)	Low dose (n=151) (%)	Inv (n=149) (%)			
Number of patients	150	150	149	147	144	146	144	144	144	138	142	141	150	(100)	150	(99)	149	(100)
Hypokalemia	7	9	17	1	0	0	0	0	0	1	0	0	9	(6)	9	(6)	17	(11)
													9		9		17	
Hypomagnesemia	3	11	3	0	0	0	0	0	0	1	0	0	2	(1)	11	(7)	3	(2)
													3		11		3	
Hypocalcemia	54	63	55	0	0	0	0	0	0	1	0	0	50	(33)	60	(40)	53	(36)
													54		63		55	
Lactic dehydrogenase increased	30	34	32	0	1	0	0	0	0	0	0	0	29	(19)	33	(22)	26	(17)
													30		35		32	
SGOT increased	43	47	42	0	1	0	0	0	0	0	0	0	39	(26)	46	(30)	40	(27)
													43		48		42	
SGPT increased	26	28	23	0	0	0	0	0	0	0	0	0	24	(16)	26	(17)	21	(14)
													26		28		23	
Other	134	119	117	0	0	0	0	1	0	0	0	0	80	(53)	82	(54)	77	(52)
													134		120		117	
Accidental injuries	1	2	4	1	0	2	0	0	1	0	0	2	2	(1)	2	(1)	6	(4)
													2		2		9	
Cardiovascular																		
Altered function	3	0	2	0	0	0	0	0	0	0	0	2	3	(2)	0	(0)	4	(3)
													3		0		4	
Altered rhythm	3	10	19	2	1	1	0	0	1	2	0	0	6	(4)	11	(7)	10	(7)
													7		11		21	
Hypercoagulability, surgical site	2	2	1	0	0	0	0	0	0	0	0	0	2	(1)	2	(1)	1	(1)
													2		2		1	
Hypercoagulability, general	3	5	2	0	0	1	0	0	0	0	0	0	3	(2)	4	(3)	3	(2)
													3		5		3	
Hypertension/hypotension	7	10	6	2	2	0	0	1	0	2	0	0	8	(5)	12	(8)	6	(4)
													11		13		6	
Other, surgical site	1	3	2	0	0	0	0	0	1	0	0	0	1	(1)	3	(2)	3	(2)
													1		3		3	
Other, general	1	1	8	0	1	1	0	0	0	0	0	0	1	(1)	2	(1)	5	(3)
													1		2		9	
Digestive system																		
Altered function	11	5	9	1	1	2	0	1	0	0	0	0	10	(7)	6	(4)	8	(5)
													12		7		11	
Bleeding episodes	0	1	0	0	0	0	0	1	0	0	0	0	0	(0)	2	(1)	0	(0)
													0		2		0	
Constipation	20	31	30	0	1	0	0	0	0	0	0	1	20	(13)	29	(19)	27	(18)
													20		32		31	
Dyspepsia	2	2	3	0	0	0	0	0	1	0	0	0	2	(1)	2	(1)	4	(3)
													2		2		4	
Irritation and Inflammation	5	4	2	2	2	3	2	0	1	0	1	0	5	(5)	7	(5)	5	(3)
													9		7		6	
Nausea and vomiting	34	30	30	5	11	4	1	2	3	2	3	4	26	(17)	30	(20)	24	(16)
													42		46		41	
Other	9	6	14	1	1	3	0	0	0	1	0	0	10	(7)	7	(5)	14	(9)
													11		7		17	
Hematic/lymphatic system																		
Anemia	85	85	86	2	1	0	0	1	1	1	0	1	79	(53)	75	(50)	74	(50)
													88		87		88	
Other	14	20	14	1	0	0	0	0	0	0	0	0	12	(8)	17	(11)	12	(8)
													15		20		14	
Immune response																		
Allergic reactions, surgical site	12	11	7	0	4	5	1	4	2	1	1	2	12	(8)	16	(11)	11	(7)
													14		20		16	
Allergic reactions, systemic	13	13	24	3	2	1	0	0	0	0	1	0	14	(9)	12	(8)	16	(11)
													16		16		25	
Autoimmune disorders	1	0	0	0	0	0	0	1	0	0	0	0	1	(1)	1	(1)	0	(0)
													1		1		0	
Infection																		
Surgical site, superficial	24	20	17	4	3	3	1	5	1	2	2	1	25	(17)	22	(15)	20	(13)
													31		30		22	
Surgical site, deep	15	10	14	2	3	0	1	2	1	0	1	1	16	(11)	13	(9)	14	(9)
													18		16		16	

Table 1: Adverse events reported for patients enrolled in clinical trial supporting approval (continued)

													# (%) of Patients Total Adverse Events		
	1st quarter postop			2nd quarter postop			3rd quarter postop			4th quarter postop					
	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl	Low dose	Inv	Ctrl (n=150) (%)	Low dose (n=151) (%)	Inv (n=149) (%)
Number of patients	150	150	149	147	144	146	144	144	144	138	142	141	150 (100)	150 (99)	149 (100)
General	8	9	18	4	0	1	2	0	1	3	1	1	13 (9)	8 (5)	15 (10)
													17	10	21
Inflammation															
Systemic	64	74	78	6	5	5	2	0	4	0	0	0	47 (31)	53 (35)	56 (38)
													72	79	87
Surgical site	72	79	67	30	24	16	6	8	10	8	3	5	78 (52)	78 (52)	71 (48)
													116	114	98
Kidney	0	1	4	0	0	0	0	0	0	0	0	0	0 (0)	1 (1)	2 (1)
													0	1	4
Liver, other	0	0	3	0	0	0	0	0	0	0	0	0	0 (0)	0 (0)	1 (1)
													0	0	3
Metabolic & nutritional	12	13	10	1	2	0	0	1	0	1	0	0	11 (7)	14 (9)	7 (5)
													14	16	10
Neuromuscular/skeletal															
Altered mental status	63	56	70	1	1	4	1	0	1	1	1	3	47 (31)	40 (26)	41 (28)
													66	58	78
Altered sensory status, general	10	8	11	0	1	3	1	0	1	0	1	0	10 (7)	5 (3)	10 (7)
													11	10	15
Altered sensory status, surgical site	16	22	27	2	4	6	2	1	0	3	1	1	18 (12)	21 (14)	28 (19)
													23	28	34
Altered motor status, general	4	7	4	1	1	0	0	0	0	0	0	0	5 (3)	8 (5)	4 (3)
													5	8	4
Altered motor status, surgical site	13	11	11	2	1	1	0	0	0	1	0	0	10 (7)	10 (7)	11 (7)
													16	12	12
Other, surgical site	160	157	158	127	92	86	33	51	37	41	33	33	121 (81)	111 (74)	109 (73)
													361	333	314
Other, general	68	59	71	9	7	19	14	5	8	11	4	7	54 (36)	43 (28)	44 (30)
													102	75	105
Respiratory system	31	45	39	0	2	1	0	1	0	0	1	0	30 (20)	38 (25)	32 (21)
													31	49	40
Skin and appendages															
Surgical area	2	1	3	1	1	1	0	0	0	0	0	0	3 (2)	2 (1)	4 (3)
													3	2	4
General	3	0	1	0	0	1	1	0	0	0	0	1	4 (3)	0 (0)	3 (2)
													4	0	3
Urogenital															
Dysfunction	6	8	10	1	0	2	1	0	1	0	1	0	7 (5)	7 (5)	8 (5)
													8	9	13
Hematuria	0	2	3	0	0	1	0	0	0	0	0	0	0 (0)	2 (1)	4 (3)
													0	2	4
Infection	1	4	4	0	0	2	0	0	1	0	0	0	1 (1)	4 (3)	6 (4)
													1	4	7

¹ No. of patients with event

² Percent of patients with event

³ Total number of events

A post-approval clinical study was designed to evaluate the performance of reamed IM nails in treating open tibia fractures in which the intramedullary canal was reamed to cortical chatter compared to the standard of care. This study was terminated prior to completion of enrollment due to safety concerns. Specifically, an increased rate of infection was observed in the rhBMP-2/ACS-treated group versus the standard of care control group (19% versus 9%, respectively). In pre-approval studies evaluating the use of rhBMP-2/ACS in combination with reamed or unreamed IM nails to treat open tibia fractures, this difference in infection rates was not observed.

Potential Adverse Events

The following is a list of potential adverse events which may occur with treatment of open tibial fractures requiring stabilization with an IM nail. Some of these adverse events may have been previously reported in the adverse events table. As with any surgery, surgical treatment of a fracture is not without risk. A variety of complications related to surgery or the use of Infuse™ Bone Graft can occur. These may occur singly or in combination. Some of these may be severe, affecting patient outcome.

- Bone fracture.
- Bowel, bladder, or gastrointestinal problems.

- Change in mental status.
- Damage to blood vessels, bleeding (which may require a blood transfusion) or cardiovascular system compromise.
- Damage to nearby tissues.
- Death.
- Development of respiratory problems.
- Disassembly, bending, breakage, loosening, and/or migration of IM nail components.
- Ectopic and/or exuberant bone formation.
- Fetal development complications.
- Foreign body (allergic) reaction.
- Incisional complications.
- Infection.
- Neurological system compromise.
- Non-union (or pseudarthrosis), delayed union, mal-union.
- Pain or discomfort.
- Rash or allergic reaction.
- Scar formation.
- Side effects from anesthesia or the surgical approach.
- Swelling.
- Tissue or nerve damage.

Note: additional surgery may be necessary to correct some of these potential adverse events.

CLINICAL RESULTS

Premarket Study Results

Safety and effectiveness of Infuse™ Bone Graft were evaluated as part of a prospective, randomized, controlled, multinational (11 countries), multi-center (49 sites) study. Subjects were randomized into one of three groups – control or one of two investigational groups (0.75 or 1.50mg/ml rhBMP-2). All subjects received wound management and fracture stabilization with an IM nail, while the investigational subjects also received Infuse™ Bone Graft at the fracture site. No restrictions were placed on the type of IM nail used or whether it was reamed or non-reamed. The use of bone wax, Gelfoam, or other collagen hemostatic agents, corticosteroids, bone growth stimulators (electrical, ultrasound, or magnetic) was specifically prohibited.

Only the subjects and an independent radiology panel were blinded with respect to treatment. The investigators were aware of the treatment assignment.

The enrolled patients had been diagnosed with an acute, open tibia fracture of Gustilo Grade I, II, IIIA or IIIB with the major component of the fracture being diaphyseal. Patients with isolated tibia fractures and those with multiple injuries were included.

Clinical and radiographic effectiveness parameters

Subjects were followed for 12 months after definitive wound closure (DWC). Evaluations were performed postoperatively at 6, 10, 14, 20, 26, 39, and 50 weeks. Adverse events, device-related or not, were evaluated over the course of the clinical trial. At each evaluation timepoint, the primary and secondary clinical and radiographic outcome parameters were evaluated. Success was determined from data collected during the initial 12 months of follow-up. Antibodies to rhBMP-2 and bovine Type I collagen were assessed preoperatively and at timepoints post-operatively. Antibodies to human Type I collagen were assessed if the antibody response to bovine Type I collagen was positive.

Clinical and radiographic fracture assessments were performed at each postoperative visit. However, the protocol did not provide the specific objective criteria that were used to determine fracture healing. The “Assessment of Fractured Limb” case report form provided for the documentation of the following parameters, but did not indicate how these were to be used to determine fracture status or how many needed to be present in order for a complete evaluation to have occurred:

- Wound (healed, not healed, not evaluated).
- Pain (absent, present, not evaluated).
- Swelling (absent, present, not evaluated).
- Tenderness (absent, present, not evaluated).
- Neurovascular status (intact, impaired, not evaluated).
- Infection (absent, present, not evaluated).
- Weight-bearing status (non-, touch down, partial, full, not evaluated).

Radiographic assessments (AP and lateral radiographs) were performed at each post-op visit. Oblique radiographs were to be used if the standard views did not adequately visualize the fracture. The radiographs were evaluated by the investigator and the radiology panel.

Investigators were to evaluate radiographic healing by assigning a score of “united”, “not united”, “uncertain union” or “uninterpretable”. No definitions for these terms were provided in the protocol.

The protocol for the independent radiology panel stated that "...fracture union was determined if there was cortical bridging and/or disappearance of the fracture lines were visible on 3 of the 4 bone aspects (anterior, posterior, medial, and lateral)...." These definitions were not available to the investigators. The first visit at which these criteria were met was considered the time of union radiographically. The independent radiographic evaluation protocol called for the review of each radiograph by all three

members of the panel. An agreement of 2 of the 3 reviewers was necessary for a determination of fracture union. The independent radiographic evaluation was performed on all available radiographs.

Adverse events were assessed for relatedness to the device and severity was based on the WHO recommendations.

Investigators were provided with the following definitions:

- "...*Nonunion* – considered to be established when a minimum of 9 months has elapsed since injury and the fracture site showed no visibly progressive signs of healing for a minimum of 3 months (no change of fracture callus)."
- *Delayed union* – insufficient fracture healing determined by radiographic and clinical assessment. A specific time point at which delayed union was defined was not provided.
- *Secondary intervention for delayed union* – any intervention, surgical or non-surgical, that was performed to induce or accelerate fracture union after DWC. Examples included use of autograft, allograft or bone graft substitutes; IM nail dynamization; exchange nailing; or noninvasive modalities, e.g., ultrasound, magnetic field, or electrical stimulation...." The decision to perform a secondary intervention for delayed union was dependent on the definition of delayed union above.

Investigators determined fracture union based on clinical judgement. The protocol did not provide the specific objective criteria that were used to determine fracture healing or deciding whether to recommend secondary interventions to promote fracture healing.

Patient demographics and accountability

The sample size estimation called for 150 subjects per treatment group. A total of 149 investigational and 150 control patients were enrolled in the study and received treatment. Only the results from the control patients and investigational patients receiving the 1.5mg/ml dose device are described below.

The demographics of the patient population were similar across all study groups except for the parameter of age, specifically the mean and range. The subjects in the investigational group were younger (mean = 33.4 years, range 18-77 years) compared to the control group (mean = 36.8 years, range = 17-87 years). Patients in the control group had a slightly larger percentage of nails that were unreamed and/or less than 9mm, while the investigational group had a higher percentage of reamed nails and /or nails that were greater than 11mm. Nail type did affect the number of secondary interventions (i.e. patients with unreamed nails had a higher incidence of secondary interventions compared to patients receiving reamed nails).

Clinical and radiographic effectiveness evaluation

The primary efficacy endpoint was defined as the proportion of subjects who required a secondary surgical intervention to promote fracture healing within 12 months of DWC. The secondary efficacy endpoints included the following:

- The proportion of subjects healed at 6 months without a secondary intervention as determined by the investigator's clinical and radiographic assessment.
- The independent radiology panel's assessment of time to fracture union.
- The pharmacoeconomic impact of the treatment.

Primary effectiveness endpoint

The rate of secondary interventions was significantly lower in the Infuse™ Bone Graft group ($p=0.001$) as described in the table below. Interventions were categorized in one of three ways – recommended by the investigator and performed, recommended by the investigator and not performed, or not recommended by the investigator but performed anyway. If any of these occurred, the patient was considered to have failed the primary endpoint and, therefore, was considered a study failure. In addition, patients who experienced screw breakage resulting in self-dynamization were also considered treatment failures.

Table 2: Number of patients with secondary interventions^a

	Control (n=150)	Investigational (n=149)
Recommended/performed	38	19
Recommended/not performed	3	5
Not recommended/performed	6	7
Self-dynamizations	19	7
Total failures/patients (%)	66/150 (44)	38/149 (26)

^aIncludes bone grafting, fibula osteotomies, exchange nailing, plate fixation, Ilizarov frame removal, external fixation placement, bone transport, IM nail dynamizations, exchange to a functional brace and electrical stimulation.

Safety and immune response evaluation

The assessment of safety consisted of an evaluation of the reported adverse events, as well as an evaluation of antibodies to rhBMP-2, bovine Type I collagen and human Type I collagen. The complete list of adverse events is described in the Adverse Events section above. Refer to this section for a description of the rates associated with infection and abnormal clinical lab values. Adverse events of special interest are discussed below.

fracture healing

The rates of hardware failure in the investigational and control groups were 18/149 (14%) and 32/150 (24%), respectively. Delayed union was the most frequent serious adverse event reported at one year; occurring in 39 (26%) control and 26 (17%) investigational patients.

The rate of nonunion was lower in the investigational group as compared to control. A total of 80/150 (53.3%) control and 56/149 (37.6%) investigational patients did not require a secondary intervention and were not radiographically healed at 12 months as determined by the independent radiology panel. For the control patients who required a secondary intervention, 18/66 (27%) reported nonunions at 12 months compared to 19/38 (50%) investigational patients. Investigational patients who required a secondary intervention were considered healed later than control patients.

abnormal bone formation

Heterotopic ossification (HO) was not a significant concern and no ectopic ossification was reported. Because only the involved tibia was evaluated radiographically, it is not clear if abnormal bone formation occurred in other anatomical locations. Although twice as many investigational patients reported hypertrophic callus formation compared to controls (8 vs. 4), no action was required to treat any of these events. A total of 12 patients experienced at least one event classified as hypertrophic callus, with the investigational group having the highest number (8) and percentage (6%) of patients with HO. No interventions were required to treat any of the HO-related events.

infections

The combined rate of deep and superficial infections of the injured limb was lower in patients with Gustilo IIIA and B fractures of the investigational group as compared to the control group [16/66 (24%) and 26/61 (43%), respectively].

immune response

The presence of antibodies was assessed using ELISA. If there was a positive response to bovine Type I collagen, the serum was also tested for antibodies to human Type I collagen. The screening ELISA cutpoint for positive antibody responses was set to 2 times the signal generated by pooled normal human sera in each ELISA. Subjects were considered to have an elevated immune response if the preoperative test was negative (titer<50) and postoperative test was positive (titer≥50) or if the preoperative test was positive and the postoperative test was positive with a three-fold higher titer than the preoperative test. There were detectable rhBMP-2 antibodies in 1 control patient and 9 investigational patients after treatment. Of the 9 investigational patients with elevated post-treatment antibody titers, 2 were elevated at visit 6 (20 weeks), the last planned assessment and data from 1 patient were unavailable for visit 6. An additional sample from 1 of these patients was collected and tested following the positive test at visit 6, and the titers decreased to <50. (No follow-up data were available for the other patient.) Anti-rhBMP-2 antibody responses were determined to be transient in 6/9 patients by 20 weeks, and in 7/9 patients after follow-up testing (samples from 2 patients were unavailable to confirm transience of the antibody response). Because of the small numbers involved, it was not possible to determine if a correlation existed between the immune response and clinical outcome.

There were 38 patients who developed antibodies to bovine Type I collagen - 9 (6%) control and 29 (20%) investigational patients. Approximately half of the patients had persistently elevated antibody titers at evaluations 20 weeks and more after DWC. Thirty categories of adverse events that may have been manifestations of an immune response were identified and all were observed to have a comparable incidence across all groups. Although there were 4 patients with an adverse event termed "allergic event", the investigators believed that there was no evidence of allergic response to the investigational treatment.

Table 3: Immune response

	Control [n (%)]	Investigational [n (%)]
Anti-rhBMP-2 antibodies	1 (1)	9 (6)
Anti-bovine Type I collagen antibodies	9 (6)	29 (20)
Anti-human Type I collagen antibodies	0 (0)	0 (0)
# healed patients with antiBMP-2 antibody response (successes)	1	6
# secondary intervention patients with antiBMP-2 antibody response (failures)	0	3

The rates of authentic antibody response to rhBMP-2 were higher than that observed for another application of rhBMP-2/ACS. When rhBMP-2/ACS was placed inside of a metallic spinal fusion cage for anterior interbody fusion treatment of degenerative disc disease, the antiBMP-2 antibody response in the investigational group was 0.7%. This compares to a 6% rate in the investigational group in the trauma study. The contribution of the trauma setting to this outcome is unknown, as is the clinical significance of the antibody response.

Post-market study results (Study 400)

As a condition of its Marketing Authorization in Europe, Wyeth Research conducted a post-market study to evaluate the safety and effectiveness of rhBMP-2/ACS (Infuse™ Bone Graft) when used specifically with reamed intramedullary (IM) nails in the treatment of open tibial shaft fractures. This study, which had a prospective, randomized, single-blind design, was conducted at 28 sites throughout Europe and South Africa. Within 14 days after the initial injury, subjects were randomized (1:1) into one of two treatment groups, which were stratified according to the Gustilo-Anderson classification of the wound at the time of the injury. All subjects received standard of care wound management and fracture stabilization with a statically-locked reamed IM nail, while the investigational subjects also received rhBMP-2/ACS (1.5 mg/ml) at the fracture site.

The study was originally designed for 300 subjects to receive one of the study treatments. A total of 277 subjects (224 men and 53 women) were randomized into one of the two treatment groups, and of these, 271 received treatment prior to the study being stopped due to an observed increase in the rate of infection in the investigational group. Study subjects ranged in age from 18 to 85 years, with an average age of 39 years. The study population was found to be comparable between treatment groups with respect to baseline demographics, morbidity, and fracture and treatment characteristics.

After receiving the assigned study treatment and being discharged from acute-care hospitalization, subjects were followed on an outpatient basis for a total of nine visits over a 52-week time period (at 6, 10, 13, 16, 20, 24, 32, 41, and 52 weeks after the baseline visit). At these protocol-specified time points, the condition of the fractured tibia and the surrounding tissues were assessed, radiographs of the fracture were evaluated, the healing status of the fracture was determined, and treatment recommendations were made. Safety evaluations were also conducted at these time points (or as reported) and included recording of adverse events (AEs) and concomitant treatments, as well as laboratory evaluations. In addition, antibody formation to BMP-2 and Type I collagen were evaluated at certain time points during the study (6, 20, and 52 weeks).

The study's primary efficacy endpoint was the proportion of subjects having a healed fracture, as demonstrated by radiographic and clinical assessment, at 13 and 20 weeks after treatment. At week 13, the rate of fracture healing was higher in the rhBMP-2/ACS group (60%) than in the SOC group (48%), but this difference failed to reach the threshold for statistical significance ($p=0.0541$). At week 20, there was no discernible difference in the proportion of subjects with healed fractures (68% in the rhBMP-2/ACS group and 67% in the SOC group).

With respect to safety, most study subjects (90%) had at least one treatment-emergent adverse event (TEAE), which is consistent with the morbidity associated with an open tibia fracture and the surgical treatment required. The most frequent TEAEs reported in the region under study were pain, peripheral edema, hardware failure, and infections. With the exception of infection, these events occurred with similar frequency in both treatment groups. Twenty-seven (27) subjects (19.4%) in the rhBMP-2 group and 15 subjects (10.9%) in the SOC group reported infections. While this difference was not statistically significant ($p=0.0645$; difference in infection risk = 0.09 [95% confidence interval, 0.0 to 0.17]), it was determined to be clinically significant by the Sponsor and was the basis for early termination of the study.

Thirty-one (22%) subjects in the rhBMP-2 group and 25 (18%) in the SOC group reported serious adverse events, four of which were deemed related to rhBMP-2. One death was reported in each treatment group during the study, but neither was thought to be related to treatment. No difference in laboratory measurements between the two groups were observed, and the rate of anti-BMP-2 and anti-bovine Type I collagen antibody formation was low and not clinically significant. No neutralizing BMP-2 antibodies and no anti-human type I collagen antibodies were found in the subjects presenting with elevated anti-BMP-2 titers and anti-bovine type I collagen antibodies, respectively.

HOW SUPPLIED

Infuse™ Bone Graft is supplied in a kit containing all the components necessary to prepare the device (i.e. the collagen sponge, a vial with the lyophilized rhBMP-2, a vial with the sterile water for reconstituting the rhBMP-2, syringes and needles).

STORAGE CONDITIONS

Store Infuse™ Bone Graft at room temperature [15–30°C (59 to 86°F)].

DOSAGE AND ADMINISTRATION

Infuse™ Bone Graft is prepared immediately prior to use from a kit containing all necessary components. Once prepared, Infuse™ Bone Graft contains rhBMP-2 at a concentration of 1.5mg/mL. The instructions for preparation must be followed and the rhBMP-2 must be reconstituted to this solution concentration of 1.5mg/ml and then distributed uniformly across the entire ACS.

Only a single Infuse™ Bone Graft kit should be used for each patient.

Infuse™ Bone Graft is implanted after the completion of IM nail fracture stabilization and wound management (i.e. at the time of soft tissue coverage). The number of Infuse™ Bone Graft kits and the volume of Infuse™ Bone Graft to be implanted is determined by the fracture anatomy. Generally, the fracture is treated with one kit. The accessible surface area of the fracture (fracture lines and defects) should be covered with Infuse™ Bone Graft. Because very few patients in the clinical trial received more than one Infuse™ Bone Graft kit (specifically, 2 kits) the response to the use of more than one kit is unknown. When determining the specific volume and placement for Infuse™ Bone Graft, the potential to induce compartment syndrome should be considered.

DIRECTIONS FOR USE

For directions for using the Infuse™ Bone Graft for tibia fractures, see the brochure entitled "Instructions for Preparation and Surgical Application."

PRODUCT COMPLAINTS

Any health care professional (e.g. customer or user of this system of products), who has any complaints or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness and/or performance of this product, should notify the distributor or Medtronic. Further, if any of the implanted INFUSE® Bone Graft ever "malfunctions," (i.e. does not meet any of its performance specifications or otherwise does not perform as intended), or is suspected of doing so, the distributor should be notified immediately. If any Medtronic product ever "malfunctions" and may have caused or contributed to the death or serious injury of a patient, the distributor should be notified immediately by telephone, fax or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report from the distributor is requested.

©2019 Medtronic Sofamor Danek USA, Inc. All rights reserved.

Covered by patents U.S. 5,013,649, U.S.5,618,924, U.S.5,166,058 and U.S.5,631,142

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©2019 Medtronic Sofamor Danek USA, Inc. Todos os direitos reservados.